

Flexible Denture Solution

High Transparency. Flexible. Pre-Polished.

High Transparency, Flexible, Lifelike Aesthetics

- Multi-shade transparent resin + high transparency printing technology create a natural look.
- Excellent rebound for a secure fit.

Prevent Deformation

Flexural Modulus

HeyGears  860 MPa

Standard Requirements¹  500 MPa

Long-lasting

Total Fracture Work

HeyGears  7000 J/m²

Standard Requirements¹  900J/m²

Tear-resistant

Elongation at break

HeyGears  50%

(No minimum requirements)²



¹ Test Method ISO 20795. ² Test Method ASTM D638

Contact the HeyGears team to learn more:

 www.heygears.com  sales@heygears.com  +1 (949) 418-9418 (USA)



Follow Us @HeyGears3D

A2D HD

Digital Workflow to Empower Your Business

- **4K high-precision** 3D printing and flexible design capabilities ensure perfect retention and comfort.
- **Extreme simplicity**, 3D printed pre-polished results, one click part removal, washing, and curing.
- Cuts hands-on time to **40** minutes per flexible partial denture (from 1-2 hours).



Great Value, Fast Delivery

- As low as **\$3** per denture base, **137** denture bases per kg of resin.
- **10 flexible dentures finished in 9 hours, from design to final product.**
- Fast turnaround, 2-day total delivery with professional design support



HeyGears Flexible Try-in

Colors

- Pink 
- Opaque Pink 
- Transparent 



Detax Denture Flex

Colors

- Pink-Transparent 
- Clear 